

Fitc Rabbit Anti-Human Tcr Γ/Δ Monoclonal Antibody - Tcr Γ/Δ

Item Number: CM0009757

Introduction

FITC-conjugated rabbit monoclonal antibody targeting human TCR γ/δ for flow cytometry. Specific for T cell receptor gamma/delta chains, ideal for immunophenotyping and T cell subset research.

[Learn More](#)

Parameter	Value
Product Name	FITC Rabbit anti-Human TCR γ/δ mAb
Remarks/Alias	TCRD; TRD@; TCRDV1; TCRG; TRG@
Species	Human
Gene ID (Human)	6964/6965
Gene ID	6964/6965
Immunogen	Recombinant protein of human TCR γ/δ
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
Shipping	Ice bag