

Fitc Rabbit Anti-Human/Monkey Cd45 Monoclonal Antibody For Fc - P08575

Item Number: CM0011945



Introduction

FITC-conjugated rabbit monoclonal antibody against human and cynomolgus CD45 (leukocyte common antigen). Validated for flow cytometry, it enables precise immune cell identification and analysis, supporting T-cell activation and signaling research.

[Learn More](#)

Parameter	Value
Product Name	FITC Rabbit anti-Human/Monkey CD45 mAb
Remarks/Alias	LCA; LY5; B220; CD45; L-CA; T200; CD45R; GP180; IMD105
Species	Human
GeneID (Human)	5788
GeneID	5788
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-577 of human CD45 (P08575)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P08575
Protein Weight	131kDa/136kDa/141kDa/143kDa/147kDa
Shipping	Ice bag