


## Fitc Rabbit Anti-Human/Monkey Cd8A Monoclonal Antibody For Fc - P01732

Item Number: CM0004144



### Introduction

FITC-conjugated rabbit monoclonal antibody specific for human and cynomolgus CD8a. Designed for flow cytometry detection of CD8+ T cells, aiding immune response and oncology research.

[Learn More](#)

| Parameter        | Value                                  |
|------------------|----------------------------------------|
| Product Name     | FITC Rabbit anti-Human/Monkey CD8a mAb |
| Remarks/Alias    | CD8; p32; Leu2; CD8alpha               |
| Species          | Human                                  |
| GeneID (Human)   | 925                                    |
| GeneID           | 925                                    |
| Immunogen        | Recombinant protein                    |
| Source           | Rabbit                                 |
| Category         | Monoclonal Antibodies                  |
| Application      | FC                                     |
| Cross Reactivity | Human, Cynomolgus                      |
| SWISS            | P01732                                 |
| Protein Weight   | 21kDa/25kDa/30kDa                      |
| Shipping         | Ice bag                                |