

Fitc Rabbit Anti-Mouse Cd32 Monoclonal Antibody For Flow Cytometry - P08101

Item Number: CM0007969

Introduction

FITC-conjugated rabbit monoclonal antibody targeting mouse CD32 (Fc gamma receptor IIB) for flow cytometry. Recognizes the extracellular domain, suitable for detecting CD32 on hemopoietic cells.

[Learn More](#)

Parameter	Value
Product Name	FITC Rabbit anti-Mouse CD32 mAb
Remarks/Alias	CD32; Fcgr2; Fcr-2; Fcr-3; Ly-17; LyM-1; Lym-1; fcRII; FcgRII; Fcgr2a; Ly-m20; Fc[g]RII; F630109E10Rik
Species	Mouse
GenelD (human)	2213
GenelD	14130
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-207 of mouse CD32 (NP_034317.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P08101
Protein Weight	31kDa/33kDa/36kDa
Shipping	Ice bag