


Fitc Rabbit Anti-Rat Cd45 Monoclonal Antibody For Fc - P04157

Item Number: CM0011923



Introduction

FITC-conjugated rabbit monoclonal antibody against rat CD45 (L-CA, T200), validated for flow cytometry. Detects leukocyte common antigen, a receptor-type tyrosine-protein phosphatase essential for T-cell activation.

[Learn More](#)

Parameter	Value
Product Name	FITC Rabbit anti-Rat CD45 mAb
Remarks/Alias	Lca; RT7; CD45; L-CA; T200
Species	Rat
GeneID (Human)	5788
GeneID	24699
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 24-546 of rat CD45 (NP_612516.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Rat
SWISS	P04157
Protein Weight	143 kDa
Shipping	Ice bag