

Fitc/5-Fam/6-Fam Rabbit Monoclonal Antibody

Item Number: CM0002056

Introduction

Rabbit monoclonal antibody targeting FITC, 5-FAM, and 6-FAM. Validated for WB and ELISA with species-independent reactivity. Ideal for detecting FITC-labeled molecules in various research and diagnostic assays.

[Learn More](#)

Parameter	Value
Product Name	FITC/5-FAM/6-FAM Rabbit mAb
Species	Human
Immunogen	Synthetic peptide, FITC
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Species independent
Shipping	Ice bag