

Fkbp1B Rabbit Pab - P68106

Item Number: CM0029527



Introduction

FKBP1B Rabbit pAb detects peptidyl-prolyl cis-trans isomerase FKBP1B (12.6 kDa FKBP) for WB, IF/ICC, ELISA in human, mouse, and rat. Suitable for immunosuppression and protein folding research.

[Learn More](#)

Parameter	Value
Product Name	FKBP1B Rabbit pAb
Remarks/Alias	OTK4; FKBP1L; PKBP1L; PPIase; FKBP12.6; FKBP1B
Species	Human
GeneID (Human)	2281
GeneID	2281
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-80 of human FKBP1B (NP_473374.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P68106
Protein Weight	12kDa
Shipping	Ice bag