

Flii Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q13045

Item Number: CM0013897

Introduction

Rabbit polyclonal antibody against human FLII (flightless-1 homolog). Validated for WB, IHC-P, ELISA; reacts with human, mouse, rat. Ideal for studies on cytoskeletal dynamics, cell adhesion, and transcription. ~145 kDa, SWISS: Q13045.

[Learn More](#)

Parameter	Specification
Product Name	FLII Rabbit pAb
Remarks/Alias	FLI; FLIL; Fli1; FLII
Species	Human
Gene ID (Human)	2314
Gene ID	2314
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 400-660 of human FLII (NP_002009.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13045
Protein Weight	145kDa
Shipping	Ice bag