

Flow Cytometry Antibody: Pe Rabbit Anti-Nhp Igd Monoclonal Antibody - Pe Rabbit Anti-Nhp Igd Mab

Item Number: CM0005839

Introduction

PE-conjugated rabbit monoclonal antibody targeting cynomolgus monkey immunoglobulin delta heavy chain. Validated for flow cytometry, cross-reacts with rhesus. Ideal for NHP IgD detection in immunology research.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-NHP IgD mAb
Species	Cynomolgus monkey
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	Flow Cytometry (FC)
Cross Reactivity	Cynomolgus, Rhesus
Shipping	Ice bag