

Fluorescein Isothiocyanate/Fitc Rabbit Pab - Anti-Fitc

Item Number: CM0021199

Introduction

Anti-FITC rabbit polyclonal antibody for WB and ELISA applications. Species-independent cross reactivity makes it ideal for universal detection of FITC-labeled molecules in IVD assay development and research.

[Learn More](#)

Parameter	Value
Product Name	Fluorescein Isothiocyanate/FITC Rabbit pAb
Immunogen	FITC
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Species independent
Shipping	Ice bag