

Fmn1 Rabbit Pab For Wb, Elisa - Q68Da7

Item Number: CM0017436

Introduction

Rabbit polyclonal antibody against human FMN1 (Formin-1) for Western blot and ELISA applications. Cross-reacts with mouse. Suitable for studying adherens junction formation and actin polymerization.

[Learn More](#)

Parameter	Value
Product Name	FMN1 Rabbit pAb
Remarks/Alias	LD; FMN; FMN1
Species	Human
Gene ID (Human)	342184
Gene ID	342184
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 350-495 of human FMN1 (NP_001096654.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q68DA7
Protein Weight	158kDa
Shipping	Ice bag