

Fndc5 Rabbit Pab - Q8Nau1

Item Number: CM0023016



Introduction

Rabbit polyclonal antibody targeting human FNDC5 (irisin). Suitable for WB, IF/ICC, and ELISA applications. Cross-reacts with human, mouse, and rat. Supports research in exercise metabolism and adipose tissue browning.

[Learn More](#)

Parameter	Value
Product Name	FNDC5 Rabbit pAb
Remarks/Alias	FRCP2; irisin; FNDC5
Species	Human
GeneID	252995
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-150 of human FNDC5 (NP_715637.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8NAU1
Protein Weight	24kDa
Shipping	Ice bag