

Formyl-Histone H2B-K120 Rabbit Monoclonal Antibody (P62807/Q16778) For Wb, Ihc-P, Elisa

Item Number: CM0000092



Introduction

Specific formyl-histone H2B-K120 rabbit monoclonal antibody for Western blot, IHC-P, and ELISA. Recognizes human, mouse, rat H2B with wide predicted cross-reactivity. Ideal for studying histone modifications and epigenetic regulation.

[Learn More](#)

Parameter	Value
Product Name	Formyl-Histone H2B-K120 Rabbit mAb
Remarks/Alias	H2B; H2BE; H2BQ; GL105; H2B.1; H2BFQ; H2BGL105; H2B-GL105; HIST2H2BE; Formyl-Histone H2B-K120
Species	Human
GeneID (Human)	8349
GeneID	3017/8349
Immunogen	A synthesized peptide derived from human Histone H2B (formyl K120)
Source	Rabbit
Category	Other Modified Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	P62807/Q16778
Protein Weight	14kDa
Shipping	Ice bag