

Foxf1 Rabbit Monoclonal Antibody For Wb And Elisa - Q12946

Item Number: CM0012022

Introduction

Reliable FOXF1 rabbit monoclonal antibody for Western blot and ELISA. Detects human and mouse FOXF1, a lung-specific transcription activator. Ideal for pulmonary development and disease research.

[Learn More](#)

Parameter	Value
Product Name	FOXF1 Rabbit mAb
Remarks/Alias	FKHL5; ACDMPV; FREAC1; FOXF1
Species	Human
GeneID (human)	2294
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human FOXF1 (Q12946)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q12946
Protein Weight	40kDa
Shipping	Ice bag