

Foxg1 Rabbit Pab - P55316

Item Number: CM0024365



Introduction

FOXG1 Rabbit pAb is a polyclonal antibody validated for Western blot and ELISA, detecting FOXG1 across human, mouse, and rat. Essential for research on brain development and transcriptional repression.

[Learn More](#)

Parameter	Value
Product Name	FOXG1 Rabbit pAb
Remarks/Alias	BF1; BF2; QIN; FKX2; HBF2; HFK1; HFK2; HFK3; KHL2; FHL3; FKHL1; FKHL2; FKHL3; FKHL4; HBF-1; HBF-2; HBF-3; FOXG1A; FOXG1B; FOXG1C; HBF-G2; FOXG1
Species	Human
GeneID (Human)	2290
GeneID	2290
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 400-489 of human FOXG1 (NP_005240.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P55316
Protein Weight	52kDa
Shipping	Ice bag