

Foxl2 Rabbit Pab For Wb, If-P, Elisa - P58012

Item Number: CM0020758



Introduction

Explore the FOXL2 Rabbit Polyclonal Antibody (PAb) validated for Western Blotting (WB), Immunofluorescence (IF-P), and ELISA. Reacts with human, mouse, rat FOXL2, a nuclear transcription factor essential for ovarian development and repressor of testis determination.

[Learn More](#)

Parameter	Value
Product Name	FOXL2 Rabbit pAb
Remarks/Alias	BPES; PFRK; POF3; BPES1; PINTO; FOXL2
Species	Human
GeneID (Human)	668
GeneID	668
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 300-376 of human FOXL2 (NP_075555.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P58012
Protein Weight	39kDa
Shipping	Ice bag