

Foxo4 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P98177

Item Number: CM0024954



Introduction

FOXO4 rabbit polyclonal antibody validated for WB, IF/ICC, ELISA with cross-reactivity to human, mouse, and rat. Ideal for investigating insulin signaling, cell cycle regulation, and transcription factor dynamics.

[Learn More](#)

Parameter	Value
Product Name	FOXO4 Rabbit pAb
Remarks/Alias	AFX; AFX1; MLLT7; FOXO4
Species	Human
GeneID (human)	4303
GeneID	4303
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P98177
Protein Weight	54kDa
Shipping	Ice bag