


Foxp1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9H334

Item Number: CM0017886



Introduction

FOXP1 rabbit polyclonal antibody suitable for Western blot, IHC-P, and ELISA. Detects endogenous FOXP1 in human, mouse, and rat. Used in transcription regulation, immunology, and stem cell research.

[Learn More](#)

Parameter	Value
Product Name	FOXP1 Rabbit pAb
Remarks/Alias	MFH; QRF1; 12CC4; hFKH1B; HSPC215; FOXP1
Species	Human
GeneID (Human)	27086
GeneID	27086
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 550-677 of human FOXP1 (NP_116071.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H334
Protein Weight	75kDa
Shipping	Ice bag