

Foxp3 Rabbit Monoclonal Antibody For Ihc-P And Elisa - Q9Bzs1

Item Number: CM0005000



Introduction

This FOXP3 rabbit monoclonal antibody is validated for IHC-P and ELISA, with cross-reactivity to human and mouse. Ideal for studying regulatory T-cell function and immune homeostasis.

[Learn More](#)

Parameter	Value
Product Name	FOXP3 Rabbit mAb
Remarks/Alias	JM2; AIID; IPEX; PIDX; XPID; DIETER
Species	Human
GeneID (Human)	50943
GeneID	50943
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human FOXP3 (NP_054728.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9BZS1
Protein Weight	47kDa
Shipping	Ice bag