

Foxp3 Rabbit Pab For Wb, Elisa - Q9Bzs1

Item Number: CM0028946



Introduction

Rabbit polyclonal antibody against human FOXP3 for WB and ELISA. Validated in human, mouse, rat. Recombinant immunogen targeting amino acids 206-300. Ideal for Treg and immune regulation research.

[Learn More](#)

Parameter	Value
Product Name	FOXP3 Rabbit pAb
Remarks/Alias	JM2; AIID; IPEX; PIDX; XPID; DIETER; FOXP3
Species	Human
GeneID	50943
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 206-300 of human FOXP3 (NP_054728.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BZS1
Protein Weight	47kDa
Shipping	Ice bag