

Fsp1/S100A4 Rabbit Pab - P26447

Item Number: CM0019867



Introduction

FSP1/S100A4 Rabbit Polyclonal Antibody validated for WB, IF/ICC, IHC-P, ELISA. Detects human and mouse S100A4, a calcium-binding protein involved in cell motility, apoptosis, and cancer metastasis research.

[Learn More](#)

Parameter	Value
Product Name	FSP1/S100A4 Rabbit pAb
Remarks/Alias	42A; 18A2; CAPL; FSP1; MTS1; P9KA; PEL98; FSP1/S100A4
Species	Human
Gene ID (Human)	6275
Gene ID	6275
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-101 of human FSP1/S100A4 (NP_002952.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P26447
Protein Weight	12 kDa
Shipping	Ice bag