

Ftcd Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - O95954

Item Number: CM0007577

Introduction

Rabbit monoclonal antibody against human FTCD, validated for WB, IF-P, IHC-P, ELISA; cross-reacts with human, mouse, rat; for studying folate metabolism and cytoskeletal organization.

[Learn More](#)

Parameter	Value
Product Name	FTCD Rabbit mAb
Remarks/Alias	LCHC1; FTCD
Species	Human
GeneID (Human)	10841
GeneID	10841
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 320-541 of human FTCD (NP_006648.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95954
Protein Weight	59kDa
Shipping	Ice bag