


Ftmt Rabbit Polyclonal Antibody (Mitochondrial Ferritin) For Wb, Elisa - Q8N4E7

Item Number: CM0016791



Introduction

High-quality rabbit polyclonal antibody against human mitochondrial ferritin (FTMT). Validated for WB and ELISA, cross-reacts with mouse. Essential for iron metabolism, mitochondrial function, and oxidative stress studies.

[Learn More](#)

Parameter	Value
Product Name	FTMT Rabbit pAb
Remarks/Alias	MTF; FTMT
Species	Human
Gene ID	94033
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 183-242 of human FTMT (NP_803431.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q8N4E7
Protein Weight	28kDa
Shipping	Ice bag