

Fto Rabbit Polyclonal Antibody - Gene Id 79068

Item Number: CM0022615

Introduction

Rabbit polyclonal antibody targeting FTO, raised against zebrafish N-terminal fragment. Validated for WB and ELISA with cross-reactivity to human, mouse, and rat. Suitable for epitranscriptomics and obesity research.

[Learn More](#)

Parameter	Value
Product Name	FTO Rabbit pAb
Species	Human
Gene ID (Human)	79068
Gene ID (Zebrafish)	553363
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-349 of zebrafish FTO (XP_009301626.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
Shipping	Ice bag