

Fuca2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Bty2

Item Number: CM0017024

Introduction

High-quality rabbit polyclonal antibody against human plasma alpha-L-fucosidase (FUCA2), validated for Western blot, IHC-P, and ELISA. Cross-reacts with human, mouse, and rat. Ideal for glycoprotein metabolism research.

[Learn More](#)

Parameter	Value
Product Name	FUCA2 Rabbit pAb
Remarks/Alias	dJ20N2.5; FUCA2
Species	Human
GeneID	2519
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 118-467 of human FUCA2 (NP_114409.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BTY2
Protein Weight	54kDa
Shipping	Ice bag