

Fuk Rabbit Polyclonal Antibody For Wb And Elisa - Q8N0W3

Item Number: CM0024549

Introduction

Rabbit polyclonal antibody against human L-fucose kinase (FUK), validated for WB and ELISA. Cross-reacts with human and mouse. Useful for fucose salvage pathway research and glycosylation disorder studies.

[Learn More](#)

Parameter	Value
Product Name	FUK Rabbit pAb
Remarks/Alias	FUK; CDGF2; 1110046B12Rik
Species	Human
Gene ID (Human)	197258
Gene ID	197258
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 795-1084 of human FUK (NP_659496.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8N0W3
Protein Weight	118kDa
Shipping	Ice bag