

Fut3 Rabbit Polyclonal Antibody (Fuct-Iii) For Wb And Elisa - P21217

Item Number: CM0016619

Introduction

Rabbit polyclonal antibody against human FUT3 (Fucosyltransferase 3), validated for WB and ELISA. Cross-reacts with rat. Ideal for Lewis blood group antigen research and glycosyltransferase studies.

[Learn More](#)

Parameter	Value
Product Name	FUT3 Rabbit pAb
Remarks/Alias	LE; Les; FT3B; CD174; FucT-III; FUT3
Species	Human
GeneID (Human)	2525
GeneID	2525
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human FUT3 (NP_000140.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P21217
Protein Weight	42kDa
Shipping	Ice bag