


Fxr1 Rabbit Pab - P51114

Item Number: CM0019683



Introduction

Rabbit polyclonal antibody targeting human FXR1 (P51114), validated for WB, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Useful for studying mRNA translation regulation, neurogenesis, and muscle development. Predicted molecular weight: 70 kDa.

[Learn More](#)

Parameter	Value
Product Name	FXR1 Rabbit pAb
Remarks/Alias	FXR1P; CMYP9A; CMYP9B; MYOPMIL; MYORIBF; FXR1
Species	Human
GeneID (Human)	8087
GeneID	8087
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 310-539 of human FXR1 (NP_001013456.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P51114
Protein Weight	70kDa
Shipping	Ice bag