

Fyn Rabbit Polyclonal Antibody For Wb, Ihc, If/Icc, Elisa - P06241

Item Number: CM0021486



Introduction

Rabbit polyclonal antibody against human Fyn (P06241), tested in WB, IHC, IF/ICC, ELISA. Recognizes human, mouse, rat Fyn. Fyn is a Src-family non-receptor tyrosine kinase involved in cell signaling, immune response, and neural function. Ideal for oncology and immunology research.

[Learn More](#)

Parameter	Value
Product Name	Fyn Rabbit pAb
Remarks/Alias	SLK; SYN; p59-FYN; Fyn
Species	Human
GeneID (human)	2534
GeneID	2534
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 200-300 of human Fyn (NP_002028.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P06241
Protein Weight	61kDa
Shipping	Ice bag