


Gaa Rabbit Monoclonal Antibody For Wb, Elisa - P10253

Item Number: CM0002218



Introduction

Anti-GAA rabbit monoclonal antibody validated for WB and ELISA. Targets human lysosomal alpha-glucosidase (P10253), essential for glycogen degradation. Ideal for Pompe disease research.

[Learn More](#)

Parameter	Value
Product Name	GAA Rabbit mAb
Remarks/Alias	LYAG; GAA
Species	Human
GeneID (Human)	2548
GeneID	2548
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 100-200 of human GAA (P10253).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P10253
Protein Weight	105kDa
Shipping	Ice bag