

Gaba Transaminase (Abat) Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P80404

Item Number: CM0004630

Introduction

Rabbit monoclonal antibody targeting human GABA transaminase (ABAT), suitable for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Recombinant immunogen (aa 1-100). Ideal for neurotransmitter degradation and mitochondrial research.

[Learn More](#)

Parameter	Value
Product Name	GABA transaminase (ABAT) Rabbit mAb
Remarks/Alias	GABAT; NPD009; GABA-AT; GABA transaminase (ABAT)
Species	Human
GeneID (human)	18
GeneID	18
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P80404
Protein Weight	56kDa
Shipping	Ice bag