


Gabbr1 Rabbit Pab - Q9Ubs5

Item Number: CM0016158



Introduction

Rabbit polyclonal antibody against human γ -aminobutyric acid type B receptor subunit 1 (GABBR1) for WB, IF/ICC, ELISA. Reacts with human, mouse, rat. Ideal for neuroscience studies on inhibitory synaptic transmission.

[Learn More](#)

Parameter	Value
Product Name	GABBR1 Rabbit pAb
Remarks/Alias	GB1; GPRC3A; GABABR1; GABBR1-3; GABBR1
Species	Human
GeneID	2550
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 200-460 of human GABBR1 (NP_068703.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UBS5
Protein Weight	108kDa
Shipping	Ice bag