


Galectin 3/Lgals3 Rabbit Monoclonal Antibody For Wb, Ihc-P, Ip, Elisa - P17931

Item Number: CM0002544



Introduction

Rabbit monoclonal antibody against human Galectin-3 (LGALS3, MAC2). Validated for WB, IHC-P, IP, ELISA. Cross-reacts with mouse and rat. Suitable for inflammation, cancer, and cell adhesion research.

[Learn More](#)

Parameter	Value
Product Name	Galectin 3/LGALS3 Rabbit mAb
Remarks/Alias	L31; GAL3; MAC2; CBP35; GALBP; GALIG; LGALS2; Galectin 3/LGALS3
Species	Human
GeneID	3958
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Galectin 3/LGALS3 (P17931).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17931
Protein Weight	26kDa
Shipping	Ice bag