

Galns Rabbit Pab For Wb - P34059

Item Number: CM0023219



Introduction

GALNS Rabbit pAb is a polyclonal antibody raised against recombinant human GALNS, validated for WB and ELISA. Detects GALNS in human, mouse, and rat samples, suitable for lysosomal storage disease research.

[Learn More](#)

Parameter	Value
Product Name	GALNS Rabbit pAb
Remarks/Alias	GAS; MPS4A; GalN6S; GALNAC6S; GALNS
Species	Human
Gene ID (Human)	2588
Gene ID	2588
Immunogen	Recombinant protein corresponding to a sequence within amino acids 27-307 of human GALNS (NP_000503.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P34059
Protein Weight	58kDa
Shipping	Ice bag