

Gamt Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q14353

Item Number: CM0030924



Introduction

GAMT Rabbit Polyclonal Antibody validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat GAMT. Targets guanidinoacetate N-methyltransferase, key for creatine synthesis. Suitable for metabolic disorder and neurodevelopmental studies.

[Learn More](#)

Parameter	Value
Product Name	GAMT Rabbit pAb
Remarks/Alias	PIG2; CCDS2; TP53I2; HEL-S-20; GAMT
Species	Human
Gene ID (Human)	2593
Gene ID	2593
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-236 of human GAMT (NP_000147.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14353
Protein Weight	26kDa
Shipping	Ice bag