

Gas2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - 043903

Item Number: CM0015327

Introduction

GAS2 Rabbit pAb detects growth arrest-specific protein 2 in WB, IHC-P, ELISA with human and mouse cross-reactivity. Recombinant immunogen yields a 35kDa band. Ideal for cytoskeleton and auditory research.

[Learn More](#)

Parameter	Value
Product Name	GAS2 Rabbit pAb
Remarks/Alias	GAS-2; GAS2
Species	Human
GeneID (Human)	2620
GeneID	2620
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human GAS2 (NP_005247.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	O43903
Protein Weight	35kDa
Shipping	Ice bag