


Gata1 Rabbit Polyclonal Antibody For Wb - P15976

Item Number: CM0014921



Introduction

Rabbit polyclonal antibody targeting human GATA1, a key erythroid transcription factor. Validated for Western blot and cross-reacts with mouse. Essential tool for research into erythropoiesis and hemoglobinopathies.

[Learn More](#)

Parameter	Value
Product Name	GATA1 Rabbit pAb
Remarks/Alias	GF1, GF-1, NFE1, XLTT, ERYF1, NF-E1, XLANP, XLTDA, GATA-1, HAEADA, GATA1
Species	Human
GeneID (Human)	2623
GeneID	2623
Immunogen	Synthetic peptide corresponding to amino acids 55-155 of human GATA1 (NP_002040.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB
Cross Reactivity	Mouse
SWISS	P15976
Protein Weight	43kDa
Shipping	Ice bag