

Gbf1 Rabbit Pab - Q92538

Item Number: CM0030798

Introduction

Rabbit polyclonal antibody against human GBF1, validated for WB, IF/ICC, ELISA. Detects human and mouse GBF1, a Golgi-specific brefeldin A-resistance guanine nucleotide exchange factor involved in COPI vesicle trafficking and Golgi assembly.

[Learn More](#)

Parameter	Value
Product Name	GBF1 Rabbit pAb
Remarks/Alias	CMT2GG; CMTDI2; CMTDIA; ARF1GEF; GBF1
Species	Human
GeneID (human)	8729
GeneID	8729
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-85 of human GBF1 (NP_004184.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q92538
Protein Weight	206kDa
Shipping	Ice bag