


Gbp1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P32455

Item Number: CM0003204



Introduction

Recombinant rabbit monoclonal antibody specific for human GBP1, validated for WB, IHC-P, ELISA. Detects human, mouse, rat GBP1. Suitable for innate immunity, autophagy, and infectious disease research.

[Learn More](#)

Parameter	Value
Product Name	GBP1 Rabbit mAb
Remarks/Alias	GBP1; guanylate-binding protein 1
Species	Human
GeneID (Human)	2633
GeneID	2633
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 493-592 of human GBP1 (P32455).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P32455
Protein Weight	68 kDa
Shipping	Ice bag