


Gbp1 Rabbit Polyclonal Antibody For Wb, Elisa - P32455

Item Number: CM0026679



Introduction

GBP1 Rabbit Polyclonal Antibody from CamelBio is designed for Western blot (WB) and ELISA applications. It targets human guanylate-binding protein 1, a key interferon-inducible GTPase involved in innate immunity and pathogen defense.

[Learn More](#)

Parameter	Value
Product Name	GBP1 Rabbit pAb
Remarks/Alias	GBP1
Species	Human
GeneID (Human)	2633
GeneID	2633
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 383-592 of human GBP1 (NP_002044.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P32455
Protein Weight	68kDa
Shipping	Ice bag