

Gcet2 Rabbit Monoclonal Antibody For Wb And Elisa - Q8N6F7

Item Number: CM0002810

Introduction

Recombinant rabbit monoclonal antibody targeting human GCET2 (HGAL, GCAT2). Validated for Western blot and ELISA, cross-reacts with human, mouse, rat. Suitable for studying germinal center B-cell biology and lymphoma.

[Learn More](#)

Core Parameter	Specification
Product Name	GCET2 Rabbit mAb
Remarks/Alias	HGAL; GCAT2; GCET2
Species	Human
Gene ID (Human)	257144
Gene ID	257144
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-178 of human GCET2 (Q8N6F7)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8N6F7
Protein Weight	21kDa
Shipping	Ice bag