


Gcn2 Rabbit Mab - Q9P2K8

Item Number: CM0007361



Introduction

GCN2 Rabbit monoclonal antibody validated for WB, IHC-P, ELISA. Reacts with human, mouse, rat. Targets eIF-2-alpha kinase GCN2, a metabolic-stress sensor regulating translation and integrated stress response.

[Learn More](#)

Parameter	Value
Product Name	GCN2 Rabbit mAb
Remarks/Alias	GCN2; PVOD2
Species	Human
GeneID (Human)	440275
GeneID	440275
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human GCN2 (NP_001013725.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9P2K8
Protein Weight	187kDa
Shipping	Ice bag