


Gemin7 Rabbit Polyclonal Antibody For Wb And Elisa - Q9H840

Item Number: CM0015950



Introduction

Rabbit polyclonal antibody targeting human GEMIN7 (SIP3), suitable for Western blot and ELISA. Useful for investigating snRNP assembly, pre-mRNA splicing, and SMN complex function. Cross-reacts with human.

[Learn More](#)

Parameter	Value
Product Name	GEMIN7 Rabbit pAb
Remarks/Alias	SIP3; GEMIN7
Species	Human
GeneID (Human)	79760
GeneID	79760
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-131 of human GEMIN7 (NP_078983.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9H840
Protein Weight	15kDa
Shipping	Ice bag