

Ggt1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P19440

Item Number: CM0022879



Introduction

Rabbit polyclonal antibody against human GGT1 (gamma-glutamyltransferase 1), validated for WB, IHC-P, and ELISA with cross-reactivity to mouse and rat. For glutathione metabolism and leukotriene research.

[Learn More](#)

Parameter	Value
Product Name	GGT1 Rabbit pAb
Remarks/Alias	GGT; GTG; GGTD; CD224; GGT 1; D22S672; D22S732; GGT1
Species	Human
Gene ID (Human)	2678
Gene ID	2678
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 370-569 of human GGT1 (NP_038265.2?)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P19440
Protein Weight	61KDa/39KDa/24KDa
Shipping	Ice bag