

GL7 Rabbit Polyclonal Antibody - Q6Yw00

Item Number: CM0021215

Introduction

Rabbit polyclonal antibody targeting Oryza sativa GL7 protein (Q6YW00, ~104 kDa). Validated for WB and ELISA, with specific reactivity in rice. Ideal for plant biology and crop research.

[Learn More](#)

Parameter	Value
Product Name	GL7 Rabbit pAb
Remarks/Alias	GL7
Species	Oryza sativa
Gene ID	4343837
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 444-570 of Oryza sativa GL7. (Q6YW00).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Oryza sativa
SWISS	Q6YW00
Protein Weight	104kDa
Shipping	Ice bag