


Glb1 Rabbit Polyclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P16278

Item Number: CM0022211



Introduction

Rabbit polyclonal antibody against human beta-galactosidase (GLB1), validated for WB, IF/ICC, IHC-P, ELISA. Cross-reacts with human, mouse, rat. Useful for lysosomal storage disease and elastogenesis research.

[Learn More](#)

Parameter	Value
Product Name	GLB1 Rabbit pAb
Remarks/Alias	EBP; ELNR1; MPS4B; GLB1
Species	Human
Gene ID (Human)	2720
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 329-677 of human GLB1 (NP_000395.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P16278
Protein Weight	76 kDa
Shipping	Ice bag