

Gli3 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P10071

Item Number: CM0009355



Introduction

High-quality Gli3 rabbit monoclonal antibody for WB, IF/ICC, ELISA. Validated in human, mouse, rat. Ideal for Hedgehog signaling and developmental biology research.

[Learn More](#)

Parameter	Value
Product Name	Gli3 Rabbit mAb
Remarks/Alias	PHS; ACLS; GCPS; PAPA; PAPB; PAP-A; PAPA1; PPDIV; GLI3FL; GLI3-190; Gli3
Species	Human
GeneID (Human)	2737
GeneID	2737
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1301-1580 of human Gli3 (P10071)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P10071
Protein Weight	170kDa
Shipping	Ice bag