


Glut3/Slc2A3 Rabbit Pab - P11169

Item Number: CM0024987



Introduction

Rabbit polyclonal antibody against GLUT3/SLC2A3, a facilitative glucose transporter, suitable for WB, IHC-P, ELISA. Cross-reacts with human, mouse, rat. Research use for glucose transport and metabolism studies. UniProt P11169.

[Learn More](#)

Parameter	Value
Product Name	GLUT3/SLC2A3 Rabbit pAb
Remarks/Alias	GLUT3; GLUT3/SLC2A3
Species	Human
Gene ID	6515
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 207-281 of human GLUT3/SLC2A3 (NP_008862.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P11169
Protein Weight	54kDa
Shipping	Ice bag