

Glutaminase (GLs) Rabbit Mab For Wb, If/Icc, If-P, Elisa - 094925

Item Number: CM0006396



Introduction

Rabbit monoclonal Glutaminase (GLS) antibody for WB, IF/ICC, IF-P, and ELISA. Reacts with human, mouse, and rat. Targets mitochondrial GLS kidney isoform, critical for glutamine metabolism. Ideal for cancer metabolism and neurotransmission research.

[Learn More](#)

Parameter	Value
Product Name	Glutaminase (GLS) Rabbit mAb
Remarks/Alias	GAC; GAM; KGA; GLS1; AAD20; DEE71; GDPAG; CASGID; EIEE71; Glutaminase (GLS)
Species	Human
GeneID (Human)	2744
GeneID	2744
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 40-140 of human Glutaminase (GLS). (NP_055720.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	094925
Protein Weight	73 kDa
Shipping	Ice bag