


Glycogen Synthase 1 (Gys1) Rabbit Monoclonal Antibody For Wb, Elisa - P13807

Item Number: CM0004965



Introduction

Rabbit monoclonal antibody against human Glycogen synthase 1 (GYS1), validated for Western blotting and ELISA. Ideal for muscle and glycogen metabolism research. Swiss-Prot: P13807.

[Learn More](#)

Parameter	Value
Product Name	Glycogen synthase 1 (GYS1) Rabbit mAb
Remarks/Alias	GSY; GYS; Glycogen synthase 1 (GYS1)
Species	Human
Gene ID (Human)	2997
Gene ID	2997
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 638-737 of human Glycogen synthase (P13807).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P13807
Protein Weight	84kDa
Shipping	Ice bag